



Technical Advancements at FRENA: India's first Nuclear Astrophysics Accelerator

FRENA: Facility for Research in Experimental Nuclear Astrophysics

PhD: University of Delhi, India Postdoc: GSI, Germany (DESPEC)

Postdoc: FRIB, USA

Akashrup Banerjee Saha Institute of Nuclear Physics, Kolkata India

Presentation Layout

• Nuclear Astrophysics: the subject

• FRENA: introduction

Some results

• Upcoming experiments

• Technical requirements and progress

• Collaboration!

Nuclear Astrophysics: the subject

- Nuclear reactions are important for energy generation in stars, stellar evolution and in nucleosynthesis.
- Most important quantity characterising these nuclear reactions: thermonuclear reaction rate at a given stellar temperature.
- This determines the rate of energy generation and the number of nuclear species created or destroyed in an astrophysical process.
- Thermonuclear reaction rate is related to the reaction cross section.
- Hence, cross section measurement is of critical importance.
- Accelerator-driven experiments try to measure reaction crosssections.




 $\frac{N_{II}}{N_{I}} = \frac{2Z_{II}}{n_{e}Z_{I}} \left(\frac{2\pi m_{e}kT}{h^{2}}\right)^{3/2} e^{-\chi/kT}$


India: accelerators

- IUAC, New Delhi: 15-16 MV terminal voltage.
- BARC-TIFR Pelletron facility: 14-15 MV terminal voltage.
- VECC: room-temperature cyclotron (K-130) and super-conducting cyclotron (K-500).

- Nuclear Astrophysics energies require lower terminal voltages.
- High-current ion beams required.
- Exceptional terminal voltage stability necessary.
- FRENA: 3 MV tandem accelerator.



Schematic Layout of FRENA






The FRENA 3 MV tandem accelerator, seen along with the ion source and high energy section






FRENA: ion sources

- Two ion injection systems- Cesium sputter source and multicusp ion injector (MCI).
- SO110: Cesium sputter source. Typically used for heavy elements.
- Programmable remotely-operable carousel with 50 cathode capacity.
- Positioned at 45° upstream of the low-energy switching magnet.
- SO120 and SO130: placed inside MCI cabinet at $\pm 30^{\circ}$. Tungsten-filament based gaseous source.
- Chopper-buncher system in light ion branch provides pulsed beams.



FRENA: voltage generation


- Parallel-fed Cockcroft-Walton principle (no moving parts).
- Oscillating signal- generated by a driver positioned outside the accelerator tank.
- Driver- transforms three-phase 440V AC supply into RF voltage.
- RF Voltage- fed into a circuit containing RF resonance coils and the dynodes which are brought into oscillation at the circuit resonance frequency.
- Oscillating potential is applied to two dynodes capacitive-coupled to the Corona rings, resulting in the generation of DC high voltage.





FRENA: charge stripping

- Electrostatic lens (Q-snout)- focuses incoming ions into the accelerator's charge stripping canal.
- Stripper canal is located in the middle of the accelerating column.
- Stripper canal is filled with Nitrogen gas.
- Nitrogen pressure of 10⁻² to 10⁻¹ mbar is maintained within the canal.
- Accelerator is filled with SF₆ gas- typical pressure of 6 bar.
- SF₆: prevents sparking and acts as a heat sink for stripper canal's turbo molecular pump.





FRENA: the accelerator

• Terminal voltage range: 200 kV to 3 MV.

• Typical beam currents

- Proton: 300 eμA

- Alpha: 70 eµA

- C: 20-50 eμA

- Pulsed beams of Proton, Alpha possible (125-4000 kHz; 1-2 ns FWHM).
- Beam energy can be varied in steps of few keV. (Useful for resonance hunting.)



FRENA: voltage stabilisation

- GVM stabilisation and Slit stabilisation modes available.
- 90° magnet (0.67 T, 1.5 m) with upstream slits to further select energy.
- 2 switching magnets (before and after 90° magnet) with 4 ports each (15°, 30°, 0°, 20°).
- Terminal voltage stability: 30 V/hour at T_{max} ; voltage ripple: ~30 V_{RMS} at T_{max} .
- Max beam power- 3 kW. Minimum transmission over entire energy range: 60%.




Akashrup Banerjee et al., Eur. Phys. J. Spec. Top. (2024) 233: 2859-2865


FRENA: terminal stability

- FRENA accelerator- tested for multiple parameters to check its suitability for nuclear astrophysics experiments.
- Stability of the terminal voltage was checked over a period of >24 hours.
- Any change in SF₆ pressure is compensated by the GVM read-back electronics-

```
GVM_{Out\, put} = Measured\, Value \times (1 + Compensation\, Factor \times SF_6\, Pressure)
```


• Compensation factor: adjustable using a 25-turn trim potentiometer.






FRENA: calibration

- A few standard experiments were carried out using DC proton beam (maximum current 100 enA).
- Neutron emission thresholds were measured.
- Gamma-resonances were studied.
- Different types of detectors were utilised: BF₃ and CeBr₃.
- Targets used: ⁷Li, ¹⁹F, ¹¹B, ²⁴Mg. Fabricated at FRENA and VECC.
- Reactions: ${}^{7}\text{Li}(p,n){}^{7}\text{Be}$, ${}^{19}\text{F}(p,n){}^{19}\text{Ne}$, ${}^{11}\text{B}(p,n){}^{11}\text{C}$, ${}^{19}\text{F}(p,\alpha\gamma){}^{16}\text{O}$, ${}^{24}\text{Mg}(p,\gamma){}^{25}\text{Al}$.
- Major collaboration between FRENA and VECC.




List of reactions used for calibration of the Terminal Voltage

Reaction studied	Resonance Energy/Reaction Threshold	Reference
¹⁹ F(p,γ) ¹⁶ O	669.0	Rev. Modern Phys. 22 (1950) 291
²⁴ Mg(p,γ) ²⁵ Al	823.0	Nucl. Phys. A 242 (1975) 519
¹⁹ F(p,γ) ¹⁶ O	872.11 0.20	Phys. Lett. 21 (1966) 61, Rev. Modern Phys. 38 (1966) 660
⁷ Li(p, n) ⁷ Be	1880.60 0.07	Phys. Lett. 21 (1966) 61, Rev. Modern Phys. 38 (1966) 660
¹¹ B(p, n) ¹¹ C	3017.17 0.07	NNDC
¹⁹ F(p, n) ¹⁹ Ne	4235.25 0.17	NNDC






Measurement of reaction threshold in ⁷Li(p, n)⁷Be




Pulse height distribution of the BF₃ detector used for the reaction threshold measurement

Total neutron cross section just above the threshold follow the variation of $(E_{set}-E_{th})^{3/2}$.

This distribution was fitted with $A(E_{set} - E_{th})$ akashrup.banerjee@saha.ac.in



To make it linear 2/3rd power of the measured neutron yield/beam charge was plotted as a function of proton energy.



FRENA: different beams

N.C. Podaru et al./Nuclear Instruments and Methods in Physics Research B 273 (2012) 231–233

Table 1 Summary of analyzed ion source beam currents and of H, He and C particle transmission through the 3 MV Tandetron™ accelerator.

Atomic species	E _{beam} (LE) (keV)	I _{beam} (analyzed) (eμΑ)	I _{beam} (Injected in Tandetron™) (eμA)	TV (kV)	Charge state (HE) (eC)	I _{beam} (HE) (eμA)	Particle transmission (%)
Н	30	648	230	200	+1	165	71.7
		469	409	2000		252	61.6
		451	421	2000		323	76.7
		620	583	3000		361	61.9
		563	524			360	68.7
		618	576			375	65.1
Не	20	63	53	2250	+2	69.5	65.6
		70	59	3000	+2	65.3	55.3
C	35	_	41.8	3000	+1	0.8	1.9
				3000	+2	15	17.9
				3000	+3	55	43.9
				3000	+4	34.6	20.7
				3000	+5	1.1	0.5
				3000	+6	0.1	0

akashrup.banerjee@saha.ac.in NISER I 12th Nov'25

233



Contents lists available at ScienceDirect

Nuclear Inst. and Methods in Physics Research, A

journal homepage: www.elsevier.com/locate/nima



The STELLA apparatus for particle-Gamma coincidence fusion

measurements with nanosecond timing

M. Heine a,b,*, S. Courtin a,b,c, G. Fruet a,b, D.G. Jenkins d, L. Morris d, D. M. Rudigier f, P. Adsley e, D. Curien a,b, S. Della Negra e, J. Lesrel e, C. Beck P. Dené a,b, F. Haas a,b, F. Hammache e, G. Heitz a,b, M. Krauth a,b, A. Meye P.H. Regan f,g, M. Richer a,b, N. de Séréville e, C. Stodel h




Fig. 4. [Color online] Photography of a 35 μ g/cm² carbon target after 50 h of exposure to 2.5 p μ A beam. The beam focus is at the upmost position so the beam spot forms a track along the outer area of the target foil close to the frame.

^a IPHC, Université de Strasbourg, Strasbourg, F-67037, France

^b CNRS, UMR7178, Strasbourg, F-67037, France

c USIAS/Université de Strasbourg, Strasbourg, F-67083, France


d University of York, York, YO105 DD, UK


e IPN d'Orsay, UMR8608, CNRS/IN2P3, PSUD 11, Orsay, F-91406, France

f Department of Physics, University of Surrey, Guildford, GU2 7XH, UK


⁸ National Physical Laboratory, Teddington, Middlesex, TW110 LW, UK


h GANIL, CEA/DSM-CNRS/IN2P3, Caen, F-14076, France






STELLA-ANDROMEDE, IJC Labs, Orsay






CASERTA, Italy



MLL, TU München

akashrup.banerjee@sa NISER I 12th Nov'25







Figure: A cartoon representation of the PERCS active-chamber. The beam and target are off-centred to protect in-situ sensors during irradiation time. The target ladder would move towards the detector cluster during counting, with the beam off.

Figure: Demonstration of the technique to be exploited with PERCS

FRENA: the possibilities

- 1-metre diameter scattering chamber commissioned, with provisions for installing detectors inside on movable arms.
- Some charged particle emitting reactions like ¹⁹F(p,α), ¹⁵N(p,α), ¹²C(¹²C,p), ¹²C(¹²C,α) etc. occurring in the CNO cycle or carbon burning phase will be measured using the scattering chamber facility.
- Charged particle detection does not suffer due to cosmic background.
- Provision for particle-gated gamma spectroscopy is possible with the chamber, with a gamma cup.



FRENA: upcoming experiments

Beam	Target	Channels of Interest
p	7Li	n,γ
p	10B	n
p	11 B	n,γ
p	12 C	γ
p	13 C	n,γ
p	14N	γ
p	19 F	η,γ,α,αγ
p	²⁷ A1	n,γ
p	⁸⁵ Rb	n
p	106,108Cd	γ
p	¹⁴⁴ Sm	γ
α	12 C	α,γ
α	16O	γ
α	¹¹³ In	α
⁷ Li	12 C	t
12C	12 C	various
16 O	16 O	various

Our Stockpile: SINP + Collaborators

Detectors

- HPGe CLOVER array (INGA)
- n-type HPGe detectors (mechanically cooled and LN₂ cooled)
- CeBr₃ detector array
- LaBr₃ detectors (different dimensions)
- NaI detector array
- Silicon detectors (various thickness)
- Scintillators- liquid, plastic
- BF₃ detectors
- ..





FRENA: developments

- Target setup suitable for high-current experiments
 - Cooled wheel for conductive targets,
 - Peltier-cooling and cryogenic-finger for radiative targets.
- Specialised target chamber for electronic veto of cosmic muons.
- Fully-automated liquid nitrogen filling system to support gamma spectroscopy campaigns at FRENA.
- Setup to support angular distribution measurement of neutron flux.



Beam Hall: 22m x 15m x 9m







Contact

frena.sinp@saha.ac.in

akashrup.banerjee@saha.ac.in

https://www.saha.ac.in/web/frenaabout-frena